

iml⁺
swiss dental implants

UNIVERSE
Catalogue



UNIVERSE
where union is strength

CONTENTS

CERTIFICATIONS

Certifications

1

UNIVERSE SYSTEM

Implant Body

3

SL surface treatment

4

Decontamination

5

Implant neck

6

Taper internal connection

6

SURGICAL KIT

CD and TD surgical Kit

8

Tools

10

GUIDED SURGERY

Guided Surgery

11

SURGICAL PLANNING

Surgical Planning

14

Preparation (Cylindrical drill surgical protocol)

16

CD Surgical Protocol

17

Preparation (Tapered drill surgical protocol)

18

TD Surgical Protocol

19

PROSTHETIC PLANNING

Lab components

21

Prosthetic parts

24

MUA

28

CAD/CAM DIGITAL DENTISTRY

CAD/CAM components

34

OVERDENTURE SOLUTIONS

FLEXATOR components

39

PACKAGING

Packaging

44

ABOUT US

About us

46

Men, materials and machines

47

Mechanical excellence

47

Quality checks

48

General sale terms

49

Scientific bibliography

50

Contacts

51

CERTIFICATIONS

Quality Maintenance is IML's guiding principle in a lifetime project embodied in the IML Production Protocol, a constantly evolving tool applied to the daily production of each component. IML SA products are marked Medical Devices:

- **FDA 510 (k) approved**

IML SA's Quality Management for the design, manufacture, and marketing of dental implants, instrumentation, and related accessories complies with the Directives and regulations in force.

- **EC (Class I) and EC 0425 (Class IIb and Class IIa),** manufactured in accordance with Medical Devices Directive 93/42/EEC and subsequent modification, amendments, and supplements.

- **[ICIM]**

UNI-EN ISO 9001:2015

UNI CEI EN ISO 13485:2012



UNIVERSE SYSTEM

Best performance in the cases:

- any bone density
- post extraction
- delayed loading
- immediate loading



Implant Body

Universe, with its particular conical shape, offers significant benefits for a wide spectrum of clinical applications providing excellent aesthetic results.

The implant features a highly innovative thread design. This special morphology, due to the combination of the spiral squared and spiral buttress threads, is able to increase the implant surface by 20% compared to the same implant without buttress threads providing a greater and more uniform area of contact between bone and implant speeding up the osteointegration process.

The alternate buttress and square double spiral loop generates a perfect balance between intrusive, compressive, and diverging forces capable of providing the bone with exceptional growth stimuli.

The strong threads are designed for the added compression and surface area required in soft bone placement in order to fill void spaces around the threads by ensuring an immediate primary stability.

Indeed, Universe, shows exceptional self-drilling properties which facilitate the implant placement and redirection in D1 and D2 bone types.


Universe stands out also for its excellent apical part design self-drilling and self-tapping.

Best performance in the cases:


- **any bone density**
- **post extraction**
- **delayed loading**
- **immediate loading**



SL surface treatment



Pic 1. Details of the cellular microstructure - branched and dendritic, with long filopodia and complex morphology - of an IML SL treated implant.



Pic 2. Details under the microscope of the surface of an IML SL treated implant.

The IML Research and Development team, in partnership with prestigious chemistry laboratories which specialise in implantology, has formulated the optimum surface treatment for their implants: the SL treatment.

The performance monitoring has been entrusted to the Polytechnic of Turin and to the University of Turin, which also periodically check production lots.

IML SL treatment is technically comparable to the best SLA® treatments, the most documented in the literature, and it is carried out using a sand-blasting technique, with different grain sizes, followed by etching of the surface using acid solutions.

The resulting surface has an appropriate structure for anchoring osteoblasts and promotes good integration of the implant with the bone tissue. In fact, this type of treatment suits any type of bone thanks to its ability to increase primary stability even in the presence of atrophic sites or compromised biological tissue.

Cell adhesion and roughness: laboratory analysis on IML SL treated implants

Aware of the key role played by the surface treatment in determining the speed and quality of osseointegration of a dental implant, IML has always invested considerable resources in designing surfaces that facilitate the cell adhesion. The tests carried out by the University of Turin on implants with IML SL treatment confirm the effectiveness of the topography

and of the surface chemistry developed by the IML Research and Development team. According to reports of the University of Turin and in accordance with protocols of the international literature, 24 hours after cell seeding on Universe SL treated implants, it is possible to appreciate the uniformity of the growth and of the cell adhesion over the entire implant surface. Furthermore, observing by microscope the nuclei (marked in blue) and the cytoskeleton (marked in red), it is evident that after 24 hours the cells not only have a very branched growth, with long filopodia and a complex morphology, but they are also multiplying in number. These are all indicators of the cellular behaviour on the IML SL treated surface. (Pic 1)

Another fundamental factor, determined by the surface treatment and constantly monitored by laboratory analysis, is the roughness, that is the result of the unevenness on the surface.

During the surface analysis on Universe implants, the CNR (National Research Council) of Turin examines the average roughness (Ra/Sa), the Skewness parameter (Rsk/Ssk), representing the prevalent symmetry, and the Kurtosis parameter (Rku/Sku), representing the indentation density. The resulting values, in relation to the international literature, confirm that the IML SL treated surfaces have an optimal roughness, homogeneously distributed. (Pic 2)

Decontamination

Even the decontamination process used for IML implants was developed in collaboration with the Research and Development team of our prestigious Italian universities partners.

This is a two-stage process, the second stage being composed of passing the implants through a plasma reactor. The "PLASMA REACTOR" project aimed to build a machine with suitable characteristics for treating dental implants and to define the optimal operating procedure and was conducted in close co-operation with the Department of Applied Science and Technology of the Polytechnic University of Turin and the Department of Surgical Sciences at the University of Turin's CIR Dental School.

Phase 1

- Objective: inorganic waste removal, mechanical machining, and surface treatments leave residues such as carbon and aluminium, universally considered possible causes of implants failing to osteointegrate;
- Procedure: liquid solution treatment;

Phase 2

- Objective: organic contamination removal, such as removal of pro-inflammatory agents;
- Procedure: treatment using gas cleaning agents applied via an electro-chemical process performed in the plasma reactor.



Implant neck



The collar, which has microgrooves, ends with a smooth switching platform that becomes more prominent as the diameter of the implant. It provides bone maintenance ensuring constant gingiva aesthetic impact.

Advantages:

- Increased bone distance in the connection between abutment and implant, the point where the bacterial load attaches
- Reduction of inflammatory phenomena
- Better peri-implant tissue preservation
- Better maintenance of crestal bone level

Taper internal connection



A true “cold welding” in the taper Implant-abutment connection is produced by contact pressure between the surface of the female cone of the implant and that of the male cone of the abutment.

The friction created between the two surfaces of equal conicity tightened to 35 Ncm generates direct, durable, and waterproof interlocking. In fact, this type of connection is the only one that approaches the ideal condition of a monoblock implant, which is the one used in the monophasic system, universally proven to be longer-lasting than biphasic systems.

However, being able to mechanically achieve a perfect Morse tapered connection requires great attention and special skills starting from the design stage.


IML's designers and mechanics know how to detect and evaluate each critical point in order to produce perfect components whose connection not only works well during the project but especially also in the patient's mouth throughout their life.

As demonstrated by numerous studies carried out around the world, the resulting system is effective and reliable.


(Bibliography page 42, rif. 1)

Advantages:


- Elimination of the passage of fluids and so bacterial colonization
- Elimination of micromovements at the interface between the components resulting in greater mechanical stability
- Low incidence of clinical complications
- Reduction of peri-implant infections
- Significant reduction in implant failures




Load 50N σ von Mises = 28.23 N/mm²



Load 150N σ von Mises = 84.86 N/mm²



Load 250N σ von Mises = 140.75 N/mm²



Load 500N σ von Mises = 273.97 N/mm²

The Universe implant has the advantage of a locking taper connection with a passing screw which achieves a precise and functionally perfect fit when coupled with a double internal hexagon. A design choice facilitating procedures during the surgical phase.

- Perfect antibacterial seal
- 12 Possible Abutment Positionings

The decisive factor in achieving complete system stability is to determine the function of the passing screw, which must guarantee the implant-abutment connection. This work is not easy as the screw is repeatedly subject to vibration fatigue which naturally tends to cause the unscrewing of the screw from its seat.

Here the solution identified by IML is the taper connection of the screw head inside the abutment.




The Universe implant connection is not only designed to guarantee correct implementation of the benefits of a locking taper connection but also those of a "one piece" implant.


In view of the studies that identify the concentration of mechanical stress in cortical areas as a possible cause of bone resorption (Bibliography page 42, rif. 2), the IML engineers have designed the Universe implant with an internal geometry that makes it possible for the implant-abutment-screw system to behave like a single one-piece implant system. Indeed the Universe implant-abutment-screw connection has a total height up to about 5.5 mm.

The benefit of this is better distribution of the load and the levers, which spread 80% of the force throughout the system instead of only centring them on the cortical area as is frequently the case in other implant systems. (Bibliography page 42, rif. 3)

SURGICAL KITS




CD surgical kit



Using a smartphone, the QR code printed on the box allows to display the surgical protocol and download it from www.iml-ch.com

TD surgical kit



There are two surgical kit of the UNIVERSE implant system: CD (cylindrical drill surgical protocol) and TD (tapered drill surgical protocol).

Surgical boxes are designed for maximum simplicity of use and made entirely of plastic materials suitable for steam sterilisation.

The instrument positions are clearly labelled in order to facilitate identification during the surgical operation and to correctly replace them after the maintenance procedure. The silicon supports secure the instruments firmly during transportation and sterilisation.

The kit contains stops that allow drills to be used safely and they are supply separately. Cylindrical drills and pilot drills are marked with indicators referring to implant height and drill stops.

All IML surgical instruments are manufactured in Surgical Steel of the highest quality that offers the best performance in terms of wear resistance and torque.

To follow carefully the directions of the surgical and prosthetic protocol and the instructions for cleaning and maintenance of the products ensures the optimal long-term performance and reliability for which products were designed.

Tools



Universe CD box for surgical instruments

	BOX-UNCD



Universe TD box for surgical instruments

	BOX-UNTD



Precision drill

		H
		U
∅	0.5	SFYS18



Cylindrical pilot drill for Universe CD

Drill ∅	2.0
	SFYS19



Cylindrical pilot drill for Universe TD

∅ Drill	2.0
	FL30-U



Drill stops kit for pilot drill (6 pcs) for Universe TD

	STF-KIT



Drill extension

	PR-FR



Guide pin

	UN-PIN



Conical drill

	Implant H					
Implant ∅	6	8	10	11.5	13	15
3.4	FC34-6	FC34-8	FC34-10	FC34-11.5	FC34-13	FC34-15
4	FC40-6	FC40-8	FC40-10	FC40-11.5	FC40-13	FC40-15
4.5	FC45-6	FC45-8	FC45-10	FC45-11.5	FC45-13	FC45-15
5	FC50-6	FC50-8	FC50-10	FC50-11.5	FC50-13	FC50-15

FC34-6 is designed for implant Ø4 h6



Red drill stops kit for drills from Ø2 to Ø4 for Universe CD

	SFYS042

Tools



Green drill stops kit for drills from Ø4.4 for Universe CD

	SFYS043

Cylindrical drill



2.7	SFYS21
3.1	SFYS22
3.6	SFYS23
3.9	SFYS20
4.4	SFYS25
4.9	SFYS55



Implant driver for motor

	H	
	10	15
	AVM-C	AVM-L



Multitool implant driver

	H		
	10	15	25
	CCIB-10	CCIB-15	CCIB-25



Dynamometric ratchet

	DN-I

Torque range: 15-45 Ncm



Fixed ratchet

	CR-U



Multitool driver for screws

	H	
	10	15
	CCIV-10	CCIV-15



Digital adapter for multitool driver

	SFYS051



Motor driver for screws

	H		
	6	12	17
	SFYS011	SFYS012	SFYS013



Smartpeg for implant

	SM-PEG



Multitool remover for abutment

	IMESTR-U1



Multitool driver for straight MUA

	AMM-U



Motor driver for straight MUA

	SFYS016



Threadformer

3.4	MC34
4.0	MC40
4.5	MC45
5.0	MC50

GUIDED SURGERY





Since 2020 IML becomes official distributor of RealGUIDE Software by 3DIEMME.

Let's see how it works:

In case a patient needs a fixed prosthesis on IML implants, his Dentist must acquire a 3D CBCT (cone beam computed tomography), a digital or conventional impressions copy and his facial esthetics parameters, sharing the data and the treatment plan idea with his Dental Technician through the cloud platform with the 3DIEMME RealGUIDE Software on Windows and Mac operating systems, as well as iPad and iPhone.

Based on the data acquired by the dentist and following the online prescription integrated in the software, the Dental Technician reconstructs a 3D model of the patient's bone and soft tissue on which the ideal virtual position of the teeth is designed. This draft project is then shared through the cloud with the implantologist who, using 3DIEMME RealGUIDE App on his iPad or iPhone, is able to proceed with the virtual positioning of dental implants, then the digital project is shared through the cloud platform with the Dental Technician who proceeds with the digital modeling of the surgical guide and the provisional prosthesis subsequently shared with the Dentist and the Implantologist. The team can discuss the clinical case through the secure chat. If the design is approved, the dental technician produces the surgical guide and the provisional prosthesis with a 3D printer and a CAD-CAM milling machine and delivers them to the implantologist who is ready to perform the surgery in a minimally invasive manner.




SURGICAL PLANNING




Surgical Planning

UNIT OF MEASUREMENT: mm


CODE	IMPLANT MEASURE (Ø x H)	A	B	C	D	E	F	G	H
		CORE Ø AT TIP	THREAD Ø AT TIP	IMPLANT Ø	INTERFACE Ø	THREAD PITCH	SURFACE TREATMENT H	SWITCHING PLATFORM H	IMPLANT H
IM34-8	3.4 X 8	1.35	2.75	3.7	3.3	1.2	7.8	0.2	8
IM34-10	3.4 X 10	1.35	2.75	3.7	3.3	1.2	9.8	0.2	10
IM34-11.5	3.4 X 11.5	1.35	2.75	3.7	3.3	1.2	11.3	0.2	11.5
IM34-13	3.4 X 13	1.35	2.75	3.7	3.3	1.2	12.8	0.2	13
IM34-15	3.4 X 15	1.35	2.75	3.7	3.3	1.2	14.8	0.2	15
IM40-6	4 X 6	1.95	2.9	4	3.6	1.2	6.1	0.9	7
IM40-8	4 X 8	1.45	3.15	4	3.6	1.2	7.8	0.2	8
IM40-10	4 X 10	1.45	3.15	4	3.6	1.2	9.8	0.2	10
IM40-11.5	4 X 11.5	1.45	3.15	4	3.6	1.2	11.3	0.2	11.5
IM40-13	4 X 13	1.45	3.15	4	3.6	1.2	12.8	0.2	13
IM40-15	4 X 15	1.45	3.15	4	3.6	1.2	14.8	0.2	15
IM45-6	4.5 X 6	1.8	3.45	4.5	3.9	1.2	6.1	0.9	7
IM45-8	4.5 X 8	1.8	3.45	4.5	3.9	1.2	7.8	0.2	8
IM45-10	4.5 X 10	1.8	3.45	4.5	3.9	1.2	9.8	0.2	10
IM45-11.5	4.5 X 11.5	1.8	3.45	4.5	3.9	1.2	11.3	0.2	11.5
IM45-13	4.5 X 13	1.8	3.45	4.5	3.9	1.2	12.8	0.2	13
IM45-15	4.5 X 15	1.8	3.45	4.5	3.9	1.2	14.8	0.2	15
IM50-6	5 X 6	2.2	3.8	5	4.2	1.3	6.1	0.9	7
IM50-8	5 X 8	2.2	3.8	5	4.2	1.3	7.8	0.2	8
IM50-10	5 X 10	2.2	3.8	5	4.2	1.3	9.8	0.2	10
IM50-11.5	5 X 11.5	2.2	3.8	5	4.2	1.3	11.3	0.2	11.5
IM50-13	5 X 13	2.2	3.8	5	4.2	1.3	12.8	0.2	13
IM50-15	5 X 15	2.2	3.8	5	4.2	1.3	14.8	0.2	15



OPTIONAL: The cover screw for bone ring can be purchased separately



Preparation (Cylindrical drill surgical protocol)



CD Surgical Protocol

The Universe CD Surgical Protocol has been developed to provide surgeons with indications on how to choose the most suitable instruments for implant site preparation, depending on the type of bone.

However, it is the duty of the surgeon to apply the most appropriate surgical protocol on the basis of his/her experience and following a thorough assessment of the clinical situation of the individual patient.

For the preparation of the implant site, IML has developed cylindrical drills with a tapered tip and depth marks in accordance with the length of the implant; they can be used with drill stops.

In case of dense D1 bone, adequate cortical bone preparation is essential in order to allow the implant to be inserted smoothly in the bone.





Fig. example of CD drilling sequence in dense bone of implant Ø3.4 h10
Consult the complete surgical protocol at www.iml-ch.com

IMPORTANT:

- Drills prepare the site 0.7 mm more than the height of the implant.
- The implant is supplied complete with cover screw
- Recommended torque max: 45 Ncm



Preparation (Tapered drill surgical protocol)



TD Surgical Protocol

The Universe TD Surgical Protocol has been developed to provide the surgeon with the most appropriate tools for bone compliance, and is also simple and practical.

The preparation of the implant site for the Universe implant is completed in 3 simple steps, after which the implant can be inserted easily:





Fig. example of TD drilling sequence of implant Ø4 h10


The preparation of the implant site is performed by tapered drills that optimise the bone available to place the implant, without waste. These drills are sized at the core of each single diameter and height of the implant to facilitate the drilling protocol reducing it to three simple steps. The particular tip shape guides the progressive advancement, respecting the bone and preparing a customized site. However, it is the duty of the surgeon to choose the most appropriate surgical protocol based on his or her experience following a thorough assessment of the individual patient's clinical situation.

IMPORTANT:

- Drills prepare the site 0.3 mm more than the height of the implant.
- The implant is supplied complete with cover screw
- Recommended torque max: 45 Ncm



1° step	2° step	3° step
<p>Precision drill</p>	<p>Pilot drill Ø2</p>	<p>Drill Ø3.4</p>



<p>Precision drill</p>	<p>Pilot drill Ø2</p>	<p>Drill Ø4.0</p>
------------------------	-----------------------	-------------------



<p>Precision drill</p>	<p>Pilot drill Ø2</p>	<p>Drill Ø4.5</p>
------------------------	-----------------------	-------------------




<p>Precision drill</p>	<p>Pilot drill Ø2</p>	<p>Drill Ø5.0</p>
------------------------	-----------------------	-------------------

PROSTHETIC PLANNING



Impression taking



Lab components



Open tray impression coping (*)

∅	U TRA-U



Slim open tray impression coping (*)

∅	U TRA-SL



Bridge open tray impression coping (*)

∅	U TRA-P



Closed tray impression coping (*)

∅	U TRA-C



Tear-off closed tray impression coping

∅	U TRA-S




Implant replica

∅	U ANI-U


















(*) Connection screw included.



Restoration



Prosthetic parts

	<p>Healing screw</p> <table border="1"> <thead> <tr> <th colspan="2"></th> <th colspan="5">H</th> </tr> <tr> <th colspan="2"></th> <th>1</th> <th>2</th> <th>3</th> <th>4</th> <th>5</th> </tr> </thead> <tbody> <tr> <td>∅</td> <td>U</td> <td>MGIU-1</td> <td>MGIU-2</td> <td>MGIU-3</td> <td>MGIU-4</td> <td>MGIU-5</td> </tr> </tbody> </table>			H							1	2	3	4	5	∅	U	MGIU-1	MGIU-2	MGIU-3	MGIU-4	MGIU-5	
		H																					
		1	2	3	4	5																	
∅	U	MGIU-1	MGIU-2	MGIU-3	MGIU-4	MGIU-5																	
	<p>Shift healing screw</p> <table border="1"> <thead> <tr> <th colspan="2"></th> <th colspan="5">H</th> </tr> <tr> <th colspan="2"></th> <th>1</th> <th>2</th> <th>3</th> <th>4</th> <th>5</th> </tr> </thead> <tbody> <tr> <td>∅</td> <td>U</td> <td>MGIU-21</td> <td>MGIU-22</td> <td>MGIU-23</td> <td>MGIU-24</td> <td>MGIU-25</td> </tr> </tbody> </table>			H							1	2	3	4	5	∅	U	MGIU-21	MGIU-22	MGIU-23	MGIU-24	MGIU-25	
		H																					
		1	2	3	4	5																	
∅	U	MGIU-21	MGIU-22	MGIU-23	MGIU-24	MGIU-25																	
	<p>Slim healing screw</p> <table border="1"> <tbody> <tr> <td>∅</td> <td>U</td> <td>MGIU-0</td> </tr> </tbody> </table>	∅	U	MGIU-0																			
∅	U	MGIU-0																					
	<p>Bridge peek healing screw (*)</p> <table border="1"> <tbody> <tr> <td>∅</td> <td>U</td> <td>MGIU-99</td> </tr> </tbody> </table>	∅	U	MGIU-99																			
∅	U	MGIU-99																					
	<p>Peek temporary straight abutment (*)</p> <table border="1"> <thead> <tr> <th colspan="2"></th> <th>H</th> </tr> <tr> <th colspan="2"></th> <th>0</th> </tr> </thead> <tbody> <tr> <td>∅</td> <td>U</td> <td>MDIU-100</td> </tr> </tbody> </table>			H			0	∅	U	MDIU-100													
		H																					
		0																					
∅	U	MDIU-100																					
	<p>15° angled peek temporary abutment (*)</p> <table border="1"> <thead> <tr> <th colspan="2"></th> <th>H</th> </tr> <tr> <th colspan="2"></th> <th>0</th> </tr> </thead> <tbody> <tr> <td>∅</td> <td>U</td> <td>MIU15-100</td> </tr> </tbody> </table>			H			0	∅	U	MIU15-100													
		H																					
		0																					
∅	U	MIU15-100																					
	<p>Temporary straight abutment (*)</p> <table border="1"> <thead> <tr> <th colspan="2"></th> <th>H</th> </tr> <tr> <th colspan="2"></th> <th>0</th> </tr> </thead> <tbody> <tr> <td>∅</td> <td>U</td> <td>MDIU-101</td> </tr> </tbody> </table>			H			0	∅	U	MDIU-101													
		H																					
		0																					
∅	U	MDIU-101																					
	<p>Spare</p> <p>Connecting screw for peek and temporary abutment</p> <table border="1"> <tbody> <tr> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td>VT-K</td> </tr> </tbody> </table>						VT-K																
		VT-K																					
	<p>Straight abutment (*)</p> <table border="1"> <thead> <tr> <th colspan="2"></th> <th>H</th> </tr> <tr> <th colspan="2"></th> <th>0</th> </tr> </thead> <tbody> <tr> <td>∅</td> <td>U</td> <td>MDIU-0</td> </tr> </tbody> </table>			H			0	∅	U	MDIU-0													
		H																					
		0																					
∅	U	MDIU-0																					

(*) Connection screw included.

Prosthetic parts



Aesthetic straight abutment (*)

		H			
		1	2	3	4
∅	U	MDIU-1	MDIU-2	MDIU-3	MDIU-4



15° angled abutment (*)

		H
		0
∅	U	MIU15-0



Shift aesthetic straight abutment (*)

		H			
		1	2	3	4
∅	U	MDIU-21	MDIU-22	MDIU-23	MDIU-24



Aesthetic 15° angled abutment (*)

		H		
		1	2	3
∅	U	MIU15-1	MIU15-2	MIU15-3



Slim straight abutment (*)

		H
		0
∅	U	MDIU-013



Aesthetic 25° angled abutment (*)

		H			
		0	1	2	3
∅	U	MIU25-0	MIU25-1	MIU25-2	MIU25-3



Flat to flat straight abutment (*)

∅	U	MDIU-99



Shift aesthetic 15° angled abutment (*)

		H		
		1	2	3
∅	U	MIU15-21	MIU15-22	MIU15-23



(*) Connection screw included.

Prosthetic parts

Non rotating Cr/Co base calcinable abutment (*)



	H	
	1	
∅	U	MDIU-70



Rotating Cr/Co base calcinable abutment (*)



	H	
	1	
∅	U	MDIU-71



Calcinable non rotating abutment (*)



	H	
	2	
∅	U	CALI-U

Connecting screw for abutment

Spare



	VT-P

(*) Connection screw included.

XL prosthetic parts

XL healing screw



	H				
	1	2	3	4	
∅	U	MGIU-XL1	MGIU-XL2	MGIU-XL3	MGIU-XL4



XL aesthetic straight abutment (*)



	H				
	1	2	3	4	
∅	U	MDIU-XL1	MDIU-XL2	MDIU-XL3	MDIU-XL4



XL aesthetic 15° angled abutment (*)




	H			
	1	2	3	
∅	U	MIU15-XL1	MIU15-XL2	MIU15-XL3











(*) Connection screw included.

MUA positioning




MUA

	<p>Straight MUA</p> <table border="1"> <thead> <tr> <th></th> <th colspan="5">H</th> </tr> <tr> <th></th> <th>1</th> <th>2</th> <th>3</th> <th>4</th> <th>5</th> </tr> </thead> <tbody> <tr> <td>∅ U</td> <td>CDIU-1</td> <td>CDIU-2</td> <td>CDIU-3</td> <td>CDIU-4</td> <td>CDIU-5</td> </tr> </tbody> </table>		H						1	2	3	4	5	∅ U	CDIU-1	CDIU-2	CDIU-3	CDIU-4	CDIU-5	
	H																			
	1	2	3	4	5															
∅ U	CDIU-1	CDIU-2	CDIU-3	CDIU-4	CDIU-5															
	<p>Shift straight MUA</p> <table border="1"> <thead> <tr> <th></th> <th colspan="5">H</th> </tr> <tr> <th></th> <th>1</th> <th>2</th> <th>3</th> <th>4</th> <th>5</th> </tr> </thead> <tbody> <tr> <td>∅ U</td> <td>CDIU-21</td> <td>CDIU-22</td> <td>CDIU-23</td> <td>CDIU-24</td> <td>CDIU-25</td> </tr> </tbody> </table>		H						1	2	3	4	5	∅ U	CDIU-21	CDIU-22	CDIU-23	CDIU-24	CDIU-25	
	H																			
	1	2	3	4	5															
∅ U	CDIU-21	CDIU-22	CDIU-23	CDIU-24	CDIU-25															
	<p>17° angled MUA (*)</p> <table border="1"> <thead> <tr> <th></th> <th colspan="2">H</th> </tr> <tr> <th></th> <th>3</th> <th>4</th> </tr> </thead> <tbody> <tr> <td>∅ U</td> <td>CIU17-3</td> <td>CIU17-4</td> </tr> </tbody> </table>		H			3	4	∅ U	CIU17-3	CIU17-4										
	H																			
	3	4																		
∅ U	CIU17-3	CIU17-4																		
	<p>30° angled MUA (*)</p> <table border="1"> <thead> <tr> <th></th> <th colspan="3">H</th> </tr> <tr> <th></th> <th>3</th> <th>4</th> <th>5</th> </tr> </thead> <tbody> <tr> <td>∅ U</td> <td>CIU30-3</td> <td>CIU30-4</td> <td>CIU30-5</td> </tr> </tbody> </table>		H				3	4	5	∅ U	CIU30-3	CIU30-4	CIU30-5							
	H																			
	3	4	5																	
∅ U	CIU30-3	CIU30-4	CIU30-5																	

(*) Connection screw included.

MUA impression taking



MUA lab components

MUA open tray impression coping (*)



	SFYP076

OPTIONAL:

Long screw for MUA impression coping



	H
	20
	SFYV011


MUA replica



	SFYP077

(*) Connection screw included.

MUA restoration




MUA prosthetic parts




MUA healing cap (*)



MUA calcinable cylinder (*)




MUA peek temporary cylinder (*)



Spare

M1.4 connecting screw for MUA prosthetic parts



Max 15 Ncm



MUA titanium cylinder (*)



MUA Cr/Co base calcinable cylinder (*)




(*) Connection screw included.

CAD-CAM DIGITAL DENTISTRY



CAD-CAM impression taking



CAD-CAM lab components



Scan body ØU (*)

Ø	U MDIU-80




CAD-CAM Implant replica ØU

Ø	U ANI-CAD

(*) Connection screw included.

CAD-CAM restoration



CAD-CAM prosthetic parts

Non rotating TBase abutment (*)

		H		
		0.5	1	2
∅	U	MDIU-50	MDIU-51	MDIU-52

Rotating TBase abutment (*)


		H		
		0.5	1	2
∅	U	MDIU-56	MDIU-53	MDIU-54

Premilled (*)

		∅	
		11.5	16
∅	U	MDIU-60	MDIU-61

(*) Connection screw included.


CAD-CAM, MUA impression taking



CAD-CAM, MUA lab components



MUA scan body (*)



CAD-CAM MUA replica (*)



Spare




CAD-CAM fixing screw for implant replica



(*) Connection screw included.

CAD-CAM, MUA restoration



CAD-CAM, MUA prosthetic parts



MUA Tbase abutment (*)

	SFYP148

Spare

M1.4 connecting screw for MUA prosthetic parts



	SFYV009


Max 15 Ncm

(*) Connection screw included.

OVERDENTURE SOLUTIONS



Overdenture solutions



Overdenture prosthetic parts, lab components



Large ball abutment

	H			
	1	2	3	4
∅ U	OTKIU-1	OTKIU-2	OTKIU-3	OTKIU-4



Flexator straight abutment

	H			
	1	2	3	4
∅ U	MDIU-201	MDIU-202	MDIU-203	MDIU-204



Titanium closed cap for large ball

	H
	4.1
	SFYP062



Flexator impression coping

	SFYP164
--	---------



Titanium open cap for large ball

	H
	2.6
	SFYP063



Flexator replica

	SFYP165
--	---------



O-ring (10 pcs)

	GOM-I
--	-------



Flexator block out spacer - white









































	SFYP162
--	---------



Flexator mid cap for lab - black

	SFYP161
--	---------

Overdenture prosthetic parts, lab components

 <p>Flexator propack 0°- 20°</p> <table border="1"><tr><td></td><td></td></tr><tr><td></td><td>SFYP166</td></tr></table>				SFYP166	 <p>Flexator mid cap HR 0°-20° - transparent</p> <table border="1"><tr><td></td><td></td></tr><tr><td></td><td>SFYP156</td></tr></table>				SFYP156
									
	SFYP166								
									
	SFYP156								
 <p>Flexator propack 20°- 40°</p> <table border="1"><tr><td></td><td></td></tr><tr><td></td><td>SFYP167</td></tr></table>				SFYP167	 <p>Flexator mid cap ZR 20°-40° - grey</p> <table border="1"><tr><td></td><td></td></tr><tr><td></td><td>SFYP157</td></tr></table>				SFYP157
									
	SFYP167								
									
	SFYP157								
 <p>Flexator titanium cap</p> <table border="1"><tr><td></td><td></td></tr><tr><td></td><td>SFYP163</td></tr></table>				SFYP163	 <p>Flexator mid cap LR 20°-40° - red</p> <table border="1"><tr><td></td><td></td></tr><tr><td></td><td>SFYP158</td></tr></table>				SFYP158
									
	SFYP163								
									
	SFYP158								
 <p>Flexator mid cap LR 0°-20° - blue</p> <table border="1"><tr><td></td><td></td></tr><tr><td></td><td>SFYP154</td></tr></table>				SFYP154	 <p>Flexator mid cap MR 20°-40° - orange</p> <table border="1"><tr><td></td><td></td></tr><tr><td></td><td>SFYP159</td></tr></table>				SFYP159
									
	SFYP154								
									
	SFYP159								
 <p>Flexator mid cap MR 0°-20° - pink</p> <table border="1"><tr><td></td><td></td></tr><tr><td></td><td>SFYP155</td></tr></table>				SFYP155	 <p>Flexator mid cap HR 20°-40° - green</p> <table border="1"><tr><td></td><td></td></tr><tr><td></td><td>SFYP160</td></tr></table>				SFYP160
									
	SFYP155								
									
	SFYP160								

Flexator tools



Flexator guide pin

	SFYS068



Flexator 3-in-1 universal driver

	SFYS067



Multitool driver for flexator
H

	6	12
	SFYS065	SFYS066



Motor driver for flexator
H

	6	12
	SFYS063	SFYS064

PACKAGING



Packaging

IML's packaging process is performed in compliance with the standards set by the EC 93/42 Directive, which guarantee the sterilisation shelf-life. The IML implants are sterilised by beta rays.

The implants are packaged in a ABS container that, in turn, is placed inside a plastic container safety seal cap. Then the plastic container is placed inside a cardboard box bearing a removable label, bearing the implant information details. Further two copies of the label are into the cardboard box, to be placed on the implant passport and on the patient's medical record sheet.




Grey ABS implant system stopper and red ABS cover screw stopper are carefully washed and dried. The dental implant is contained in titanium spacers.



The transparent grey fumè Polypropylene (PP) container is closed with a white Polypropylene (PP) stopper with a safety seal.



The cardboard box (3.5 x 6.2 x 3.5 cm) must be stored in a dry place at room temperature.



ABOUT US

The Swiss company IML SA Swiss Dental Implants was founded in 2009 by a close-knit team of professionals with twenty-year experience in the dental industry, especially dealing with implants.

Its engineers continuously strive to find effective solutions for new implantology needs, ones that meet the expectations of the most demanding professionals.

Main aim: to offer oral implantology that is Simple, Safe and Stable through time.

These "3Ss" summarise the guidelines the Company has established for its own standards and are pursued in every action it takes every day.



Men, materials and machines

Only the best raw materials, the most advanced technology, and the best professional.

These secrets of IML guarantee excellent products, free from manufacturing defects.

- Super-skilled operators able to develop a man-machine relationship able to optimise the features of their tools to achieve maximum performance
- Top quality titanium for medical use. grade 4 for implants and grade 5 for prostheses. IML titanium is exclusively imported from the United States, is guaranteed free from manufacturing defects and radioactivity
- Mechanical production using latest generation sliding head machines

Mechanical excellence

How important is it for the mechanical work in the connection of an implant or in the head of a screw to be well-executed?

Just as important as it is that the abutment remains well screwed to the implant.

IML is fully aware of the issues generated by all types of production defects and knows how to resolve them, and above all, it knows how to obtain, and systematically repeat, a PERFECT MECHANICAL EXECUTION.

For example, IML guarantees 5 thousandths of a millimeter tolerance on the measure of the hexagonal connection of the implant, on every single implant.



Quality checks

Control of quality or quality control? A play on words, useful in explaining that checking is not enough for IML.


Control in IML is synonymous with uncompromising elimination of all those components that present the slightest imperfection even if only aesthetic.

It means making a commitment to selling only very specific components in order to be "as precise as the Swiss".

It means that we must fully take on the cost of this commitment both pursued and maintained ethically and proudly by IML and by taking the patient's health and the surgeon's skills into consideration.

Process:

- Identification of each individual component's critical points
- Drafting documents with a list of the critical points specific to each individual component indicating the sequence of checks to be carried out
- Over 30 checks are performed on 100% of the components manufactured in the various manufacturing phases:
 - Dimensional controls
 - Removal of burrs and dross
 - Functional tests to remove non-perfect components are performed on 100% of the components
- The operator signs off each check to certify that he or she accepts responsibility for the checks made
- Regular laboratory analyses check conformity of implant surfaces



General sale terms

VALIDITY

This catalogue is the 2019_1 edition and replaces any previous editions.

EXCLUSION OF RESPONSIBILITY

The dental implants manufactured by IML SA (hereinafter also "IML") and other IML medical devices may be used only with original IML components and instruments, following the instructions provided inside the package. The use of devices manufactured by other companies or manufactured by IML, but not belonging to the same implant line shall invalidate the warranty and terminate any explicit or implicit obligation of IML.

The clinical protocols provide the practitioner with a reference guide and shall not be construed as an alternative to the user's training and professional experience.

The practitioner using IML medical devices must ensure that device being used is suited to the patient and to existing circumstances. IML does not acknowledge any explicit or implicit responsibility, nor shall it have any responsibility for any direct, indirect, criminal or other damage deriving from or associated with any mistakes in professional judgements, in the practical application or insertion of IML products. The practitioner is also under the obligation to keep up-to-date at all times with the most recent developments and applications of IML medical devices.

Any descriptions of the products, depictions, illustrations in catalogues, illustrative and sales materials, price lists or other informative documents issued and distributed by IML in any form are provided only for explicative purposes. The Purchaser acknowledges that he/she does not purchase the products on the basis of such descriptions and/or illustrations.

PRODUCTS MODIFICATIONS

The images of the products shown in the catalogues, and in all IML publications, are for illustrative purposes only.

IML reserves the right, at any time and at its sole discretion, to make changes to the products, codes and descriptions without obligation of notice.

COPYRIGHT AND TRADEMARKS

Protection of privacy pursuant to law. It is strictly prohibited to copy or publish all or any part of this catalogue without the written authorisation of IML in any paper and electronic form.

ORDER PLACEMENT

IML dental implants and medical devices can be ordered by quoting the product code and desired quantity, in the following ways:

Telephone: 0041 (0) 916001310 Mon – Fri 09:00-13:00 / 14:00-17:00, E-mail: info@immediateload.com

PRODUCT SHIPMENTS

Except in the event of force majeure, goods will be shipped to the address indicated by the client according to the terms indicated by the supplier at the time of acceptance of the order.

RETURNED MATERIALS

The right to have purchased products replaced may be exercised within 5 working days from the products delivery date and must be pre-authorized by IML in writing. The right to replacement shall be cancelled in the event the product integrity (packaging and/or its content) is compromised. For example in the cases in which IML observes:

- that the package has been opened
- that the package has been damaged (even if still sealed)
- that the product has been damaged other than for transportation.

If IML does not accept to replace the product, it will return it, charging shipment costs to the Purchaser.

PRICES

The values indicated in the official price list are net of VAT, transport costs and any bank charges for payment collection. Unless previously agreed otherwise, prices charged will be the ones indicated in the official price list in force at the time of the order acceptance.

CHANGES TO THE PRODUCTS AND PRICES

IML reserves the right to make changes to the products contained in this catalogue and to their prices or to discontinue their production at any time, without having to give prior notice.

APPLICABLE LAW AND JURISDICTION

These Sales Terms shall be governed and interpreted according to Swiss Civil Code and the exclusive place of jurisdiction shall be the courts of Lugano.

Scientific bibliography

Connessione Conometrica Morse-taper connection

- Di Carlo F, Marincola M, Quaranta A, Bedini R, Pecci R (2008) Analisi MicroTac di impianti a connessione conometrica. *DentCadmos*76(3):1-6.
- Bedini R, Ioppolo P, Pecci R, Rizzo F, Di Carlo F, Quaranta M (2007) Studio in vitro sulla connessione di sistemi implantari dentali. *Rapporti ISTISAN* 07/7.
- Hermann JS, Schoolfield JD, Schenk RK, Buser D, Cochran DL (2001) Influence of the size of the microgap on crestal bone changes around titanium implants. A histometric evaluation of unloaded non-submerged implants in the canine mandible. *J Periodontol* 72:1372-1383.
- Crespi R, Cappare P, Gherlone E. Radiographic evaluation of marginal bone levels around platform-switched and non-platform-switched implants used in an immediate loading protocol. *Int J Oral Maxillofac Implants* 2009.
- Dibart S, Warbington M, Su MF, Skobe Z (2005) In vitro evaluation of the implant-abutment bacterial seal: the locking taper system. *Int J Oral Maxillofac Implants* 20:732-737.
- Persson LG, Leckholm U, Leonhardt A, Dahlen G, Lindhe J (1996) Bacterial colonization on internal surfaces of Branemark system implant components. *Clin Oral Implants Res* 7:90-95.
- Hansson S. Implant-Abutment Interface: Biomechanical study of Flat Top versus Conical. 2000: Biomechanics, department of polymeric materials, Chalmers University of Technology Goteborg, Sweden.
- Dibart S. et al. In vitro evaluation for the implant-abutment bacterial seal: The locking Taper system. *Int. J. Oral Maxillo-facial Implants.* 2005; 20: 732-737.

Riassorbimento osseo Bone reabsorbition

- Horowitz, R., Current Implant Designs to Maintain Crestal Bone and Gingiva, *Functional Esthetics & Restorative Dentistry: Series 1, Number 2, Dental Implants*, p. 88-90, 2008.
- Kitamura E, Stegaroiu R, Nomura S, Miyakawa O. Influence of marginal bone resorption on stress around an implant—A three-dimensional finite element analysis. *J Oral Rehabil* 2005; 32:279–286.
- King GN, Hermann JS, Schoolfield JD, Buser D, Cochran DL. Influence of the size of the microgap on crestal bone levels in non-submerged dental implants: a radiographic study in the canine mandible. *J Periodontol.* 2002;73(10):1111-7.
- Tawil G, Aboujaoude N, Youman R. Influence of prosthetic parameters on the survival and complication rates of short implants. *Int J Oral Maxillofac implants* 2006; 21(2): 275-282. 5) Leonard, G., Coelho, P., Polyzois, I., Stassen, L., Claffey, N., A study of the bone healing kinetics of plateau versus screw root design *titanium dental implants. Clinical Oral Implants Research*, 2009, 20 (3), 232-239.
- Tada S, Stegaroiu R, Kitamura E, Miyakawa O, Kusakari H. Influence of implant design and bone quality on stress/strain distribution in bone around implants: A 3-dimensional finite element analysis. *Int J Oral Maxillofac Implants* 2003; 18:357–368.
- Bozkaya, D., and Müftü, S., Muftu, A., Evaluation of Load Transfer Characteristics of Five Different Implant Systems in Compact Bone at Different Load Levels by Finite Element Analysis. *Journal of Prosthetic Dentistry*, Vol. 92 No.6, p.523-530, December 2004.
- Bidez MW, Misch CE. Issues in bone mechanics related to oral implants. *Implant Dent* 1992; 1(4): 289-294. 9) Petrie CS, Williams JL. Shape optimization of dental implant designs under oblique loading using the p-version finite element method. *J Mechanics Med Biol* 2002;2: 339–345.
- Lemons, J.E., Biomaterials, Biomechanics, Tissue Healing, and Immediate-Function Dental Implants, *Journal of Oral Implantology*, Vol. XXX No. 5, 2004.
- Venuleo, C., Chuang, S.K., Weed, M., Dibart, S., Long term bone level stability on Short Implants: A radiographic follow up study. *Indian Journal of Maxillofacial and Oral Surgery*, September 2008, Vol. 7: No.3, p. 340-345.
- Mericske-Stern R, Grutter L, Rosch R, et al: Clinical evaluation and prosthetic complications of single tooth replacements by non-submerged implants. *Clinical Oral Implants Res* 2001; 12: 309-318.
- Dibart S, Warbington M, Fan Su M, et al: In vitro evaluation of the implant-abutment bacterial seal: the Locking Taper system. *Int J Oral Maxillofacial Implants* 2005; 20: 732-737.
- Li Shi, Alex S. L. Fox. Shape Optimization of Dental Implants. *Int J Oral Maxillofac Implants.* 22:911-920; 2007. 15) Renouard F, Nisand D, Impact of implant length and diameter on survival rates. *Clinical Oral Implants Res.* 2006 Oct;17 Suppl 2:35-51.
- Schulte, J., Flores, A., Weed, M., Crown-to-implant ratios of single tooth implant-supported restorations, *Journal of Prosthetic Dentistry*, Vol 98, Issue 1, July 2007, Pages 1 UNI EN 1642. Dispositivi medici per l'odontoiatria. Impianti dentali. Milano: Ente Nazionale Italiano di Unificazione; 1997. ISO TR 11175. Dental implants. Guidelines for development dental implants. Geneva: International Organization for Standardization; 1993.
- Brunski JB. Biomaterials and biomechanics in dental implant design. *Int J Oral Maxillofac Implants* 1988; 3(2): 85-97.
- Lekholm U, Zarb GA. Tissue integrated prostheses: osseointegration in clinical dentistry. Chicago: Brånemark, Zarb & Albrektsson Eds.; 1985.
- Albrektsson T, Brånemark PI, Hansson HA, Lindström J. Osseointegrated titanium implants. Requirements for ensuring a long lasting, direct bone anchorage in man. *Acta Orthop Scand* 1981; 52(2): 155-70.
- Adell R, Lekholm U, Rockler B, Brånemark PI. A 15-year study of osseointegrated implants in the treatment of the edentulous jaw. *Int J Oral Surg* 1981; 10(6): 387-416.
- Brånemark PI, Adell R, Breine U, Hansson BO, Lindström J, Ohlsson A. Intra-osseous anchorage of dental prostheses I. Experimental studies. *Scand J Plast Reconstr Surg* 1969; 3: 81-100.
- Asikainen P, Klemetti E, Vuillemin T, Sutter F, Rainio V, Kotilainen R. Titanium implants and lateral forces. An experimental study with sheep. *Clin Oral Implants Res* 1997; 8(6): 465-468.
- Keltjens HMAM, Creugers TJ, Creugers NHJ. Three different filling materials in overdenture abutments; a 30-months evaluation. *J Dent Res* 1997; 76(5):1103.
- Weinaender M. Bone growth around dental implants. *Dent Clin North Am* 1991; 35: 585-601
- M. Marincola, L. Paracchini, V. Morgan, J. Schulte. Impianti corti: principi biomeccanici e predittività a lungo termine. *Quintessenza Internazionale* 2008, Settembre-Ottobre 2008, 45-53, 24, 5bis.
- M. Danza, I. Zollino, L. Paracchini, G. Riccardo, S. Fanali, F. Carinci. 3D finite element analysis to detect stress distribution: Spiral family implants. *J. Maxillofac. Oral Surg* 8(4): 334-339.

Carichi e leve Loading and leverage

- Brunski JB. Biomaterials and biomechanics in dental implant design. *Int J Oral Maxillofac Implants.* 1988;3(2):85–97.
- Albrektsson T, Brånemark PI, Hansson HA, Lindström J. Osseointegrated titanium implants. Requirements for ensuring a long lasting, direct bone anchorage in man. *Acta Orthop Scand.* 1981;52(2):155–70.
- Adell R, Lekholm U, Rockler B, Brånemark PI. A 15-year study of osseointegrated implants in the treatment of the edentulous jaw. *Int J Oral Surg.* 1981;10(6):387–416.
- Brånemark PI, Adell R, Breine U, Hansson BO, Lindström J, Ohlsson A. Intra-osseous anchorage of dental prostheses I. Experimental studies. *Scand J Plast Reconstr Surg.* 1969;3:81–100.
- Asikainen P, Klemetti E, Vuillemin T, Sutter F, Rainio V, Kotilainen R. Titanium implants and lateral forces. An experimental study with sheep. *Clin Oral Implants Res.* 1997;8(6):465–468.
- Keltjens HMAM, Creugers TJ, Creugers NHJ. Three different filling materials in overdenture abutments; a 30-months evaluation. *J Dent Res.* 1997;76(5):1103.
- Danza M, Zollino I, Paracchini L, Riccardo G, Fanali S, Carinci F. 3D finite element analysis to detect stress distribution: Spiral family implants. *J Maxillofac Oral Surg.* 8(4):334–339.
- Danza M, Zollino I, Paracchini L, Voza I, Guidi R, Carinci F. 3D finite element analysis comparing standard and reverse conical neck implants: Bone platform switching. *EDI Journal.* 2010;2:334–339.
- Danza M, Quaranta A, Carinci F, Paracchini L, Pompa G, Voza I. Biomechanical evaluation of dental implants in D1 and D4 bone by Finite Element Analysis. *Minerva Stomatol.* 2010;59(6):305–13.
- Danza M, Paracchini L, Carinci F. Analisi agli elementi finiti per la definizione della distribuzione degli stress meccanici negli impianti. *Dental Cadmos.* 2012;80(10):598–602.
- Marrelli M, Maletta C, Inchingolo F, Alfano M, Tatullo M. Three-point bending tests of zirconia core/veneer ceramics for dental restorations. *Int J Dent.* 2013;2013:831976. doi:10.1155/2013/831976.
- Inchingolo F, Tatullo M, Abenavoli FM, Marrelli M, Inchingolo AD, Palladino A, Inchingolo AM, Dipalma G. Oral piercing and oral diseases: a long term retrospective study. *Int J Med Sci.* 2011;8(8):649–52.
- Paduano F, Marrelli M, White LJ, Shakesheff KM, Tatullo M. Odontogenic Differentiation of Human Dental Pulp Stem Cells on Hydrogel Scaffolds Derived from Decellularized Bone Extracellular Matrix and Collagen Type I. *PLoS One.* 2016 Feb 16;11(2):e0148225. doi:10.1371/journal.pone.0148225.
- Marrelli M, Falisi G, Apicella A, Apicella D, Amantea M, Cielo A, Bonanome L, Palmieri F, Santacroce L, Giannini S, Di Fabrizio E, Rastelli C, Gargari M, Cuda G, Paduano F, Tatullo M. Behaviour of dental pulp stem cells on different types of innovative mesoporous and nanoporous silicon scaffolds with different functionalizations of the surfaces. *J Biol Regul Homeost Agents.* 2015 Oct–Dec;29(4):991–7.
- Tatullo M, Marrelli M, Falisi G, Rastelli C, Palmieri F, Gargari M, Zavan B, Paduano F, Benagiano V. Mechanical influence of tissue culture plates and extracellular matrix on mesenchymal stem cell behavior: A topical review. *Int J Immunopathol Pharmacol.* 2016 Mar;29(1):3–8. doi:10.1177/0394632015617951 Review.
- Inchingolo F, Tatullo M, Marrelli M, Inchingolo AM, Inchingolo AD, Dipalma G, Flace P, Girolamo F, Tarullo A, Laino L, Sabatini R, Abbinante A, Cagiano R. Regenerative surgery performed with platelet-rich plasma used in sinus lift elevation before dental implant surgery: an useful aid in healing and regeneration of bone tissue. *Eur Rev Med Pharmacol Sci.* 2012 Sep;16(9):1222–6.
- Marrelli M, Tatullo M. Influence of PRF in the healing of bone and gingival tissues. *Clinical and histological evaluations.* *Eur Rev Med Pharmacol Sci.* 2013 Jul;17(14):1958–62.
- Inchingolo F, Tatullo M, Pacifici A, Gargari M, Inchingolo AD, Inchingolo AM, Dipalma G, Marrelli M, Abenavoli FM, Pacifici L. Use of dermal-fat grafts in the post-oncological reconstructive surgery of atrophies in the zygomatic region: clinical evaluations in the patients undergone to previous radiation therapy. *Head Face Med.* 2012 Dec 5;8:33. doi:10.1186/1746-160X-8-33. Review. • Gargari M, Prete V, Pujia A, Ceruso FM. Full-arch maxillary rehabilitation fixed on 6 implants. *Oral Implantol (Rome).* 2013 Jul 15;6(1):1–4.
- Carinci F, Brunelli G, Franco M, Viscioni A, Rigo L, Guidi R, Strohmenger L. A retrospective study on 287 implants installed in resorbed maxillae grafted with fresh frozen allogeneous bone. *Clin Implant Dent Relat Res.* 2010 Jun 1;12(2):91–8.
- Scarano A, Murmura G, Carinci F, Lauritano D. Immediately loaded small-diameter dental implants: Evaluation of retention, stability and comfort for the edentulous patient. *European Journal of Inflammation.* 2012;10(1):19–23.
- Fanali S, Carinci F, Zollino I, Brugnati C, Lauritano D. One-piece implants installed in restored mandible: A retrospective study. *European Journal of Inflammation.* 2012;10(1):37–41.
- Andreasi Bassi M, Andrisani C, Lico S, Ormanier Z, Ottria L, Gargari M. Guided bone regeneration via a preformed titanium foil: Clinical, histological and histomorphometric outcome of a case series. *Oral Implantol (Rome).* 2016 Oct–Dec;9(4):164–174.



Contacts

IML SA

Administrative and operational location:
Via Moree, 16 - Ingresso B
6850 Mendrisio (Switzerland)
Tel: +41 (0)916001310
www.iml-ch.com
info@immediateload.com

COME AND VISIT US

You are most welcome to come visit us on a guided tour of our Company.
Do not hesitate to contact us for a date.

Distance by car from airports:

- Lugano Agno (LUG) - Switzerland > 22 km - 25 min
- Milan Malpensa (MXP) - Italy > 57 km - 45 min
- Milan Linate (LIN) - Italy > 72 km - 60 min
- Milan Orio al Serio (BGY) Italy > 102 km - 90 min

imt⁺
swiss dental implants

